

## IFNgamma Gene Promoter Reporter Vector

### Product Insert

**Cat. No. LR1011**


**Shelf Life** One year from date of receipt under proper storage conditions.

**Storage** Store at -20°C

### Contents

Gene Promoter Reporter Vector	20ul of 500ng/ul; 10ug
Empty Control Reporter Vector	20ul of 500ng/ul; 10ug
Sequence information is available on our website	

Gene Promoter	IFN gamma
Description	Designed to study regulation of promoter for Homo sapiens interferon, gamma (IFNG)
RefSeq ID (for downstream gene)	NM_000619



**Sequence (5'-3')**

```

CCCTGGGAATATTCTCTACACTGTATTTCAAGGATTTAATATGACAAAAAGAATGTCAAATACCTTATTAACAATGTAGT
ATATTGATGCATACTGAAGTACTATTTGGGATATATTGGTTAAATACAATATATTTAAATATATTTACCTTTTAAA
AAAACTTTTATTAATGAGGCTACTAGATCATTAAATTTACCTGTGTGGCTTGTATTGTATTCTACTGGGCAGTGCTGA
TCTAGAGCAATTTGAAACTTGTGGTAGATATTTACTAACCAACTCTGATGAAGGACTTCCTCACCAAAATTTGCTTTTTA
ACCGCATTCTTTCCCTTGCTTTCTGGTCATTTGCAAGAAAAATTTAAAAGGCTGCCCTTTGTAAAGGTTTGAGAGGCC
TAGAATTTCTGTTTTCACTTGTTCACCAACCACAAGCAAATGATCAATGTGCTTTGTGAATGAAGAGTCAACATTTACCA
GGCGAAAGTGGGGAGGTACAAAAAATTTCCAGTCCTTGAATGGTGTGAAGTAAAAGTGCCTTCAAAGAATCCCACCAGA
ATGGCACAGGTGGGCATAATGGGTCTGTCTCATCGTCAAAGGACCCAAGGAGTCTAAAGGAACTCTAACTACAACACCC
AAATGCCACAAAACCTTAGTTATTAATACAACTATCATCCCTGCCTATCTGTCACCATCTCATCTTAAAAAaCTTGTGA
AAATACGTAATCCTCAGGAGACTTCAATTAGGTATAAATACCAGCAGCCAGAGGAGGTGCAGCACATTGTTCTGATCATC
TGAAGATCAGCTATTAGAAGAGAAAGATCAGTTAAGTCCTTTGGACCTGATCAGCTTGATACAAGAACTACTGATTTCAA
CTTCTTTGGCTTAATTCTCTCGGAAAC
  
```

---

**Contacting  
Panomics**

For ordering information or technical support, contact the appropriate resource provided below according to your geographical location.

<b>Location</b>	<b>U.S. Corporate Headquarters</b>
Address	Panomics, Inc. 6519 Dumbarton Circle Fremont CA 94555 USA
Telephone	1.510.818.2600
FAX	1.510.818.2610
Email	info@panomics.com
Technical Support	<a href="mailto:techsupport@panomics.com">techsupport@panomics.com</a> or <a href="tel:1.877.726.6642">1.877.726.6642 option 3</a>
Ordering Information	orders@panomics.com

<b>Location</b>	<b>European Headquarters</b>
Address	Panomics Srl Via Sardegna 1 20060 Vignate Milano Italy
Telephone	+39.02.95.360.250
FAX	+39.02.95.360.992
Email	info_europe@panomics.com
Technical Support	<a href="mailto:techsupport_europe@panomics.com">techsupport_europe@panomics.com</a>
Ordering Information	order_europe@panomics.com

<b>Location</b>	<b>Asia Pacific Headquarters</b>
Address	Panomics, Inc. 16F Gemdale Plaza Tower A, No. 91 Jianguo Road Beijing 100022 P.R. China
Telephone	+86.10.59208157
FAX	+86.10.59208111
Email	info_asia@panomics.com
Technical Support	<a href="mailto:techsupport_asia@panomics.com">techsupport_asia@panomics.com</a>
Ordering Information	order_asia@panomics.com

For an updated list of FAQs and product support literature, visit our website at [www.panomics.com](http://www.panomics.com).

---